

VITAL®

Canon
CANON GROUP

Non-DICOM Conformance Statement EasyViz 8.3

September 22. 2021
Revision 1.3
Document ID: 2020.02.001

Vital Images A/S
Krumtappen 4, Etage 3
DK-2500 Valby
Denmark

www.vitalimages.com/contact

Vital Images shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Vital Images reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Vital Images products.

Prepared By: Vital Images

Table of Contents

1	Non-DICOM Conformance Statement Overview	3
2	Specification	3
3	Supported MIME types	3
4	Annotations on non-DICOM images	4
5	No export or transfer of MHD data.....	4

Revision History

<i>Revision</i>	<i>Date</i>	<i>Author</i>	<i>Changes</i>
1.0	March 20, 2020	Karen Piech	New for 8.0
1.1	September 25, 2020	Karen Piech	Updated for 8.1
1.2	March 24, 2021	Karen Piech	Updated for 8.2
1.3	September 22, 2021	Karen Piech	Updated for 8.3

1 Non-DICOM Conformance Statement Overview

The EasyViz image display implements the necessary IHE Mobile access to Health Documents (MHD) services necessary to query, retrieve and display non-DICOM images and other text-based documents.

2 Specification

Non-DICOM content is obtained in accordance with the Find Document References [ITI-67] specification using the following parameters.

<i>ITI-67 parameter name</i>
patient or patient.identifier

3 Supported MIME types

Support for non-DICOM content includes both image and text-based MIME types.

<i>Extension</i>	<i>Description</i>	<i>MIME type</i>
.bmp	Windows OS/2 Bitmap Graphics	image/bmp image/x-windows-bmp
.gif	Graphics Interchange Format (GIF)	image/gif
.htm .html	HyperText Markup Language (HTML)	txt/html
.ico	Icon format	image/x-icon
.jpeg .jpg	JPEG images	image/jpeg image/jpg image/pjpeg
.pbm	Portable BitMap	image/x-portable-bitmap
.pdf	Portable Document Format (PDF)	application/pdf
.pgm	Portable GrayMap	image/x-portable-graymap image/x-portable-greymap
.png	Portable Network Graphics	image/png
.pnm	Portable Any Map	image/x-portable-anymap

<i>Extension</i>	<i>Description</i>	<i>MIME type</i>
.ppm	Portable Pixmap	Image/x-portable-pixmap
.tif .tiff	Tagged Image File Format (TIFF)	image/tiff image/x-tiff
.txt	Text	text/plain
.wpg	WordPerfect Graphics (WPG) file	application/wpg application/x-wpg
.xbm	X BitMap	image/x-xbitmap image/xbm image/x-xbm
.xpm	X Pixmap	image/x-ppixmap image/xpm image/x-xpm

4 Annotations on non-DICOM images

EasyViz supports the addition of annotations on non-DICOM images. Measurements are calculated in pixels, however, and they cannot be saved with the image, even as a compound presentation.

5 No export or transfer of MHD data

EasyViz currently does not support the export or transfer of MHD data:

- It is not possible to export MHD from the CAM to a CD-ROM or DVD.
- It is not possible to transfer MHD from the Workflow Manager or from the Patient Study Browser.