

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Informatics	Vitrea Connection	9.1	12/4/23

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Antepartum Education	Content Consumer	View Option
Antepartum Education	Content Consumer	None
Antepartum History and Physical	Content Consumer	View Option
Antepartum History and Physical	Content Consumer	None
Antepartum Laboratory	Content Consumer	View Option
Antepartum Laboratory	Content Consumer	None
Antepartum Summary	Content Consumer	View Option
Antepartum Summary	Content Consumer	None
Anatomic Pathology Structured Reports	Content Consumer	View Option
Anatomic Pathology Structured Reports	Content Consumer	None
Pathology Workflow	Image Manager / Image Archive	None
Access to Radiology Information	Image Manager / Image Archive	None
Audit Trail and Node Authentication	Audit Record Repository	None
Audit Trail and Node Authentication	Audit Record Repository	ATX: TLS Syslog

Internet address for vendor's IHE information :
<https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/>

Links to Standards Conformance Statements for the Implementation

HL7	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/
DICOM	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Informatics	Vitrea Connection	9.1	12/4/23

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Audit Trail and Node Authentication	Secure Application	Radiology Audit Trail
Audit Trail and Node Authentication	Secure Application	None
Audit Trail and Node Authentication	Secure Application	ATX: TLS Syslog
Audit Trail and Node Authentication	Secure Application	STX:TLS 1.2 Floor using BCP195
Audit Trail and Node Authentication	Secure Node	Radiology Audit Trail
Audit Trail and Node Authentication	Secure Node	None
Audit Trail and Node Authentication	Secure Node	ATX: TLS Syslog
Audit Trail and Node Authentication	Secure Node	STX:TLS 1.2 Floor using BCP195
Basic Image Review	Image Manager / Image Archive	None
Basic Patient Privacy Consents	Content Consumer	Basic Patient Privacy Acknowledgement View
Basic Patient Privacy Consents	Content Consumer	None
Basic Patient Privacy Consents	Content Creator	None
Basic Patient Privacy Consents	Content Creator	Basic Patient Privacy Acknowledgment
Cardiac Imaging Report Content	Content Consumer	View Option

Internet address for vendor's IHE information :
<https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/>

Links to Standards Conformance Statements for the Implementation

HL7	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/
DICOM	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Informatics	Vitrea Connection	9.1	12/4/23

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cardiac Imaging Report Content	Content Consumer	None
Consistent Presentation of Images	Image Manager / Image Archive	None
Cath Report Content	Content Consumer	View Option
Cath Report Content	Content Consumer	None
Consistent Time	Time Client	Secure NTP.
Consistent Time	Time Client	None
Composite Triage and Nursing Note	Content Consumer	View Option
Composite Triage and Nursing Note	Content Consumer	None
Chest X-Ray CAD Display	Image Manager / Image Archive	None
Digital Breast Tomosynthesis	Image Manager / Image Archive	None
MR Diffusion Imaging	Image Manager / Image Archive	None
Displayable Reports	Enterprise Report Repository	None
Encounter-Based Imaging Workflow	Encounter Manager	None
Evidence Documents	Image Manager / Image Archive	None
ED Physician Note	Content Consumer	View Option

Internet address for vendor's IHE information :
<https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/>

Links to Standards Conformance Statements for the Implementation

HL7	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/
DICOM	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Informatics	Vitrea Connection	9.1	12/4/23

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
ED Physician Note	Content Consumer	None
Emergency Department Referral	Content Consumer	View Option
Emergency Department Referral	Content Consumer	None
E-Nursing Summary	Content Consumer	View Option
E-Nursing Summary	Content Consumer	None
EMS Transport Summary	Content Consumer	View Option
EMS Transport Summary	Content Consumer	None
Image Fusion	Image Manager / Image Archive	None
Immunization Content	Content Consumer	View Option
Immunization Content	Content Consumer	None
Imaging Object Change Management	Change Requester	None
Imaging Object Change Management	Image Manager / Image Archive	None
Interfacility Transport Summary	Content Consumer	View Option
Interfacility Transport Summary	Content Consumer	None
Key Image Note	Image Manager / Image Archive	None

Internet address for vendor's IHE information :
<https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/>

Links to Standards Conformance Statements for the Implementation

HL7	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/
DICOM	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Informatics	Vitrea Connection	9.1	12/4/23

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Labor and Delivery History and Physical	Content Consumer	View Option
Labor and Delivery History and Physical	Content Consumer	None
Labor and Delivery Summary	Content Consumer	View Option
Labor and Delivery Summary	Content Consumer	None
Mammography Image	Image Manager / Image Archive	None
Maternal Discharge Summary	Content Consumer	View Option
Maternal Discharge Summary	Content Consumer	None
Multi-Patient Queries	Document Consumer	None
Multi-Patient Queries	Document Registry	None
Newborn Discharge Summary	Content Consumer	View Option
Newborn Discharge Summary	Content Consumer	None
NM Image	Image Manager / Image Archive	None
Nursing Note	Content Consumer	View Option
Nursing Note	Content Consumer	None

Internet address for vendor's IHE information :
<https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/>

Links to Standards Conformance Statements for the Implementation

HL7	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/
DICOM	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Informatics	Vitrea Connection	9.1	12/4/23

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Patient Demographics Query	Patient Demographics Consumer	None
Patient Demographics Query	Patient Demographics Supplier	None
CT/MR Perfusion with Imaging Contrast	Image Manager / Image Archive	None
Patient Information Reconciliation	Image Manager / Image Archive	HL7 v2.5.1
Patient Information Reconciliation	Image Manager / Image Archive	None
Patient Information Reconciliation	Performed Procedure Step Manager	None
Patient Identifier Cross-referencing	Patient Identity Source	None
Patient Identifier Cross-referencing	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-referencing	Patient Identifier Cross-reference Manager	None
Patient Plan of Care	Content Consumer	View Option
Patient Plan of Care	Content Consumer	None
Postpartum Visit Summary	Content Consumer	View Option
Postpartum Visit Summary	Content Consumer	None

Internet address for vendor's IHE information :
<https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/>

Links to Standards Conformance Statements for the Implementation

HL7	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/
DICOM	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Informatics	Vitrea Connection	9.1	12/4/23

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Patient Care Plan	Content Consumer	View Option
Patient Care Plan	Content Consumer	None
Simple Image and Numeric Report	Enterprise Report Repository	None
Scheduled Workflow	Image Manager / Image Archive	HL7 v2.5.1
Scheduled Workflow	Image Manager / Image Archive	None
Scheduled Workflow	Performed Procedure Step Manager	None
Scheduled Workflow.b	Image Manager / Image Archive	None
Scheduled Workflow.b	Performed Procedure Step Manager	None
Triage Note	Content Consumer	View Option
Triage Note	Content Consumer	None
Sharing Lab Report	Content Consumer	View Option
Sharing Lab Report	Content Consumer	None
Cross-Enterprise Document Reliable Interchange	Document Recipient	None
Cross-Enterprise Document Reliable Interchange	Document Source	None

Internet address for vendor's IHE information :
<https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/>

Links to Standards Conformance Statements for the Implementation

HL7	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/
DICOM	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Informatics	Vitrea Connection	9.1	12/4/23

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross-Enterprise Document Sharing	Document Administrator	DEPRECATED 2017: Update Metadata
Cross-Enterprise Document Sharing	Document Administrator	None
Cross-Enterprise Document Sharing	Document Administrator	DEPRECATED:Delete Metadata
Cross-Enterprise Document Sharing	Document Consumer	Basic Patient Privacy Proof
Cross-Enterprise Document Sharing	Document Consumer	None
Cross-Enterprise Document Sharing	Document Consumer	Basic Patient Privacy Enforcement
Cross-Enterprise Document Sharing	Document Registry	Document Metadata Update
Cross-Enterprise Document Sharing	Document Registry	None
Cross-Enterprise Document Sharing	Document Registry	Patient Identity Feed (HL7 V2)
Cross-Enterprise Document Sharing	Document Source	Folder Management
Cross-Enterprise Document Sharing	Document Source	None
Cross-Enterprise Document Sharing	Document Source	Document Addendum
Cross-Enterprise Document Sharing	Document Source	Document Transformation
Cross-Enterprise Document Sharing	Document Source	Document Replacement

Internet address for vendor's IHE information :
<https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/>

Links to Standards Conformance Statements for the Implementation

HL7	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/
DICOM	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/

Link to general information on IHE :www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Informatics	Vitrea Connection	9.1	12/4/23

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross-Enterprise Document Sharing	Document Source	Basic Patient Privacy Enforcement
Cross-Enterprise Document Sharing for Imaging	Document Consumer	None
Cross-Enterprise Document Sharing for Imaging	Document Registry	None
Cross-Enterprise Document Sharing for Imaging	Document Repository	None
Cross-Enterprise Document Sharing for Imaging	Imaging Document Consumer	None
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	Set of DICOM Instances
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	None
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	PDF Report
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	CDA Wrapped Text Report

Internet address for vendor's IHE information :
<https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/>

Links to Standards Conformance Statements for the Implementation

HL7	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/
DICOM	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Informatics	Vitrea Connection	9.1	12/4/23

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	CDA Imaging Report with Structured Headings
Cross-Enterprise Sharing of Medical Summaries	Content Consumer	None
Cross-Enterprise Sharing of Medical Summaries	Content Consumer	View Option
Cross-Enterprise Sharing of Medical Summaries	Content Creator	Discharge Summary
Cross-Enterprise Sharing of Medical Summaries	Content Creator	None
Cross-Enterprise Document Sharing of Scanned Documents	Content Consumer	View Option
Cross-Enterprise Document Sharing of Scanned Documents	Content Consumer	None
Cross-Enterprise Document Workflow	Content Consumer	None
Cross-Enterprise Document Workflow	Content Consumer	View Option
Cross-Enterprise Document Workflow	Content Creator	None

Internet address for vendor's IHE information :
<https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/>

Links to Standards Conformance Statements for the Implementation

HL7	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/
DICOM	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Informatics	Vitrea Connection	9.1	12/4/23

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross-Enterprise Document Workflow	Content Updater	View Option
Cross-Enterprise Document Workflow	Content Updater	None
Exchange of Personal Health Record Content	Content Consumer	View Option
Exchange of Personal Health Record Content	Content Consumer	None
XAD-PID Change Management	Document Registry	None
XAD-PID Change Management	Patient Identifier Cross-reference Manager	None

Internet address for vendor's IHE information :
<https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/>

Links to Standards Conformance Statements for the Implementation

HL7	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/
DICOM	https://www.mi.medical.canon/support-and-education/customer-support/ethics-and-compliance/compliance-documents/

Link to general information on IHE : www.ihe.net