

| <b>IHE Integration Statement</b>  |   |  |             |
|---|---|--|-------------|
| <b>Vendor</b>   | <b>Product Name</b>   | <b>Version</b>   | <b>Date</b> |
| Vital Images, Inc   | Vitrea  | V7.11  | 16-May-19   |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:                        |   |  |             |
| <b>Integration Profiles Implemented</b>   | <b>Actors Implemented</b>   | <b>Options Implemented</b>   |             |
| <b>Radiology</b> IHE Radiology Technical Framework (Rev 17.0)   |   |  |             |
| Access to Radiology Information   | Image Display   | Multiple Sources   |             |
| Scheduled Workflow  | Image Display   | None   |             |
| Portable Data for Imaging   | Image Display   | None   |             |
|   | Portable Media Creator  | None   |             |
|   | Portable Media Importer   | None   |             |
| <b>IT Infrastructure</b> IHE IT Infrastructure Technical Framework (Rev 15.0)   |   |  |             |
| Consistent Time   | Time Client   | None   |             |
| Internet address for vendor's IHE information:<br><a href="http://www.vitalimages.com/about-us/compliance-documents">http://www.vitalimages.com/about-us/compliance-documents</a> |   |  |             |
| <b>Link to Standards Conformance Statements for the Implementation</b>  |   |  |             |
| <b>DICOM</b>  | <a href="http://www.vitalimages.com/about-us/compliance-documents">http://www.vitalimages.com/about-us/compliance-documents</a> |  |             |
| <b>Links to General Information on IHE</b>  |   |  |             |
| In North America:<br><a href="http://www.ihe.net/">http://www.ihe.net/</a>  | In Europe:<br><a href="http://www.ihe.net/">http://www.ihe.net/</a>   | In Japan:<br><a href="http://www.ihe-i.org/">http://www.ihe-i.org/</a> |             |