

IHE Integration Statement			
Vendor	Product Name	Version	Date
Vital Images, Inc	Vitrea	V7.14	08-Aug-20
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology IHE Radiology Technical Framework (Rev 17.0)			
Access to Radiology Information	Image Display	Multiple Sources	
Scheduled Workflow	Image Display	None	
Portable Data for Imaging	Image Display	None	
	Portable Media Creator	None	
	Portable Media Importer	None	
IT Infrastructure IHE IT Infrastructure Technical Framework (Rev 15.0)			
Consistent Time	Time Client	None	
Internet address for vendor's IHE information: http://www.vitalimages.com/about-us/compliance-documents			
Link to Standards Conformance Statements for the Implementation			
DICOM	http://www.vitalimages.com/about-us/compliance-documents		
Links to General Information on IHE			
In North America: http://www.ihe.net/	In Europe: http://www.ihe.net/	In Japan: http://www.ihe-j.org/	