

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Canon Medical Informatics, Inc	Vitrea	V7.15	21-Apr-22
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Radiology</b> IHE Radiology Technical Framework (Rev 17.0)			
Access to Radiology Information	Image Display	Multiple Sources	
Scheduled Workflow	Image Display	None	
Portable Data for Imaging	Image Display	None	
	Portable Media Creator	None	
	Portable Media Importer	None	
<b>IT Infrastructure</b> IHE IT Infrastructure Technical Framework (Rev 15.0)			
Consistent Time	Time Client	None	
Internet address for vendor's IHE information: <a href="http://www.vitalimages.com/about-us/compliance-documents">http://www.vitalimages.com/about-us/compliance-documents</a>			
<b>Link to Standards Conformance Statements for the Implementation</b>			
<b>DICOM</b>	<a href="http://www.vitalimages.com/about-us/compliance-documents">http://www.vitalimages.com/about-us/compliance-documents</a>		
<b>Links to General Information on IHE</b>			
In North America: <a href="http://www.ihe.net/">http://www.ihe.net/</a>	In Europe: <a href="http://www.ihe.net/">http://www.ihe.net/</a>	In Japan: <a href="http://www.ihe-j.org/">http://www.ihe-j.org/</a>	